

Indications,Safety,and Warnings

Disposal bipolar forceps

Indications

Disposable bipolar forceps are used for diagnostic and therapeutic in clinical laparoscopic and endoscopic surgeries and must not be used for other purposes. They are equipped with a high-frequency generator, enabling grasping, coagulation of human tissues.

Potential complication

The most common potential complication is thermal injury, and other complications include, but are not limited to Fever, Bleeding, hematoma, Pain, death

Warnings and precautions

Disposable bipolar forceps is provided sterile and is intended for use in a single procedure only.

Discard after use. Do not resterilize.

Please refer to instructions for use for complete contraindication and risk information.